

EN

YLED-1F @ Acrobat® Swing

MAVIG

Lamp Housing
with integrated
wide-range power supply,
protection class IP44

Sterilisable Handle
for convenient positioning
and focusing



Model: LE5117100 YLED-1F, focusable LED lamp with sterilisable handle, integrated power supply and Acrobat® Swing ceiling column, including Acrobat® Swing extension/spring arm

Made in Germany

CE, UL

YLED-1F – LED Lamp with Acrobat® Swing Ceiling System

The YLED-1F covers a wide range of applications, from diagnostic to medical examination procedures. It is intuitively controlled by means of the side-mounted control panel and the sterilisable handle.

The compact YLED-1F contains 17 powerful LED modules. All modules provide a combined light intensity of 70,000 Lux. In addition, a wide range power supply is integrated into the housing.

- **Central light intensity** (at 1 m / 39.4" distance) / **Colour temperature: 70,000 lx / 4100 ± 200 K**
- **Colour rendering index R_a at 4100 Kelvin (CRI): 95**
- **Focusable light field size** (at 1 m / 39.4" distance): **140 - 250 mm / 5.5 - 9.8"**
- **Electronic brightness control: 50 - 100 %** (incl. endoscopy light function)
- **Sterilisable handle**
- **Temperature increase in head area: 0.5 K**
- **Power consumption (total): 24 VA, at mains voltage 100 - 240 VAC, 50 - 60 Hz**
- **Long lifetime of LEDs: ≥ 60.000 h**
- **Working area: 70 - 140 cm / 27.6 - 55.1" – Height adjustment** (on Acrobat Swing spring arm) **aprox. 123 cm / 48.4"**
- **Housing colour / Weight YLED-1F: RAL 9002 / 2.4 kg**
- **Hazardous substances (EU Directive 2011/65/EU) – RoHS compliant**
- **Housing - Protected against splashed water IP44 / Fire protection class V0**
- **Medical Products Directive 93/42/EEC / Use according to DIN VDE 0100-710 / Approvals: CE, UL**

MAVIG GmbH

X-Ray Protection and Medical
Suspension Systems
Stahlgruberring 5
81829 Munich
Germany
PO Box 82 03 62
81803 Munich
Germany
Phone +49 (0) 89 / 420 96 - 0
Fax: + 49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com
www.mavig.com

ТОВ „ОНКОМЕДИКА“

**Повноважений
представник в Україні**
01014, Україна, м. Київ, вул.
Звіринська, 63
+38 (044) 228 84 17
+38 (097) 959 10 89
info@oncomedica.ua
www.oncomedica.ua
www.mavig.com.ua