


EN

RP278

MAVIG



Models:	RP2781L	Ovarian Protection, Large	RP2781G	Ovarian Protection Set,
	RP2781M	Ovarian Protection, Medium		all three sizes
	RP2781S	Ovarian Protection, Small		

Made in Germany

CE 0302

RP278 – Ovarian X-Ray Protection for Female Supine Patients

The concept of this MAVIG ovarian shield is unsurpassed and meets all the requirements of radiological practice. Radiation exposure of the female gonads is significantly reduced without compromising the diagnostic technique.

For optimal fine adjustment, each model can be varied again by means of the modular principle. The RP278 is, of course, easy to clean and disinfect. A complete set* offers you the best possible protection for all your supine patients.

- **Protective Material** Standard lead
- **Lead Equivalent** 1.00 mm Pb
- **Outer Cover/Colour** Soft PVC cover in light grey
- **Sizes**
 - Large – Height 12 cm
 - Medium – Height 9 cm
 - Small – Height 7 cm

* Also available as the complete set (all 3 sizes) with price saving = RP2781G

MAVIG GmbH

X-Ray Protection and Medical
Suspension Systems
Stahlgruberring 5
81829 Munich
Germany
PO Box 82 03 62
81803 Munich
Germany
Phone +49 (0) 89 / 420 96 - 0
Fax: + 49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com
www.mavig.com

ТОВ „ОНКОМЕДІКА“
Уповноважений
представник в Україні
01014, Україна, м. Київ, вул.
Звіринецька, 63
+38 (044) 228 84 17
+38 (097) 796 71 30
+38 (097) 959 10 89
info@oncomedica.ua
www.oncomedica.ua
www.mavig.com.ua