

EN

AW426


MAVIG

Optional to AW426

Basket AW426K for accessories

**Optional to AW426**

Holder AW426S for accessories



AW42610 for 10 pieces of clothing

Article-No.:	AW42605	Mobile storage system 5 pieces of clothing
	AW42610	Mobile storage system 10 pieces of clothing

AW426 – Mobile System for X-Ray Protective Clothing

The AW426 mobile system is easy, comfortable and safe to move. The combination of two hanger parts allow costumes in particular to be hung up quickly, easily and stored properly. The offset of the hangers means that coats and front aprons can also be stored carefully.

- **Functional metal construction**
- **For X-ray protection clothing** (5 or 10 pieces of clothing)
- **Optional: Accessories basket** (Article-No. AW426K, max. two baskets can be attached on each side)
Accessories holder (Article-No. AW426S, one holder can be attached on each side)
- **Suitable and gentle storage**
- **Easy assembly**
- **Technical Data**

Weight (net)	44.0 kg (AW42605), 62.0 kg (AW42610)
Mobility	4 heavy load casters, 2 with brakes
Storage capacity	5 or 10 X-ray protective aprons/costumes/coats (max. 75.0 kg for AW42605 / 150.0 kg for AW42610)
Outer dimensions	AW42605: 1550 x 595 x 970 mm (H x W x D) AW42610: 1550 x 1100 x 970 mm (H x W x D) AW426K: 70 x 305 x 240 mm (H x W x D)
Colour	Chrom / RAL 9010

MAVIG GmbH

X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich Germany

Phone +49 (0) 89 / 420 96 - 0 Fax
+49 (0) 89 / 420 96 - 200 e-Mail
info@mavig.com

**ТОВ „ОНКОМЕДИКА“
Уповноважений представник
в Україні**

01014, Україна, м. Київ, вул.
Звіринецька, 63
+ 38 (044) 228 84 17
+ 38 (097) 959 10 89
info@oncomedica.ua
www.oncomedica.ua